
Product Datasheet

Rabbit Anti-GLS1 antibody orb157166-CF647

Description	Rabbit polyclonal antibody to GLS1
Reactivity	Bovine, Canine, Equine, Human, Mouse, Porcine, Rabbit, Rat, Sheep
Conjugation	Unconjugated
Custom conjugations	Unconjugated, Biotin, FITC, HRP, CF647, CF488A, CF594, CF405M
Tested Applications	ELISA, ICC, IF, IHC, WB
Storage	Shipped at 4°. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Expiration Date	12 months from date of receipt.