

Product Datasheet

Glutaredoxin 5 antibody {CF647} orb157192-CF647

Description	Rabbit polyclonal antibody to Glutaredoxin 5.
Reactivity	Human, Mouse, Rat
Predicted Reactivity	Equine,Bovine,Porcine,Sheep
Conjugation	CF647
Tested Applications	ICC, IF
Concentration	100 µg in 500 µl (0.20 mg/ml)
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Dilution Range	To be determined by the end user
Expiration Date	12 months from date of receipt.