

Product Datasheet

KMT6 (Phospho-Ser21) antibody {CF647} orb157761-CF647

Description	Rabbit polyclonal antibody to KMT6 (Phospho-Ser21).
Reactivity	Human, Mouse, Rat
Predicted Reactivity	Canine, Equine, Bovine, Porcine, Sheep, Zebrafish, Xenopus
Conjugation	CF647
Tested Applications	ICC, IF
Concentration	100 µg in 500 µl (0.20 mg/ml)
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Dilution Range	To be determined by the end user
Expiration Date	12 months from date of receipt.