

Product Datasheet

Rabbit Anti-Mitoferrin 2/SLC25A28 antibody orb157881-CF647

Description	Rabbit polyclonal antibody to Mitoferrin 2/SLC25A28
Reactivity	Bovine, Canine, Mouse, Porcine, Rat, Sheep
Conjugation	Unconjugated
Custom conjugations	Unconjugated, Biotin, FITC, HRP, CF647, CF488A, CF594, CF405M
Tested Applications	ELISA, ICC, IF, IHC, WB
Storage	Shipped at 4°. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Expiration Date	12 months from date of receipt.