

## Product Datasheet

### NUMBL antibody {CF647} orb158040-CF647

|                             |                                                                                  |
|-----------------------------|----------------------------------------------------------------------------------|
| <b>Description</b>          | Rabbit polyclonal antibody to NUMBL.                                             |
| <b>Reactivity</b>           | Human, Mouse, Rat                                                                |
| <b>Predicted Reactivity</b> | Canine, Equine, Bovine, Sheep                                                    |
| <b>Conjugation</b>          | CF647                                                                            |
| <b>Tested Applications</b>  | ICC, IF                                                                          |
| <b>Concentration</b>        | 100 µg in 500 µl (0.20 mg/ml)                                                    |
| <b>Storage</b>              | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| <b>Note</b>                 | For research use only.                                                           |
| <b>Dilution Range</b>       | To be determined by the end user                                                 |
| <b>Expiration Date</b>      | 12 months from date of receipt.                                                  |