

Product Datasheet

P2RX5 antibody {CF647} orb158076-CF647

Description	Rabbit polyclonal antibody to P2RX5.
Reactivity	Human, Mouse, Rat
Predicted Reactivity	Canine,Bovine,Porcine,Sheep
Conjugation	CF647
Tested Applications	ICC, IF
Concentration	100 µg in 500 µl (0.20 mg/ml)
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Dilution Range	To be determined by the end user
Expiration Date	12 months from date of receipt.