

Product Datasheet

Rabbit Anti-STRA6 antibody orb158514-CF647

Description	Rabbit polyclonal antibody to STRA6
Reactivity	Canine, Equine, Mouse, Rabbit, Rat
Conjugation	Unconjugated
Custom conjugations	Unconjugated, Biotin, FITC, HRP, CF647, CF488A, CF594, CF405M
Tested Applications	ELISA, ICC, IF, IHC
Storage	Shipped at 4°. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Expiration Date	12 months from date of receipt.