

Product Datasheet

TTC21B antibody {CF647} orb158662-CF647

Description	Rabbit polyclonal antibody to TTC21B.
Reactivity	Human, Mouse, Rat
Predicted Reactivity	Canine, Equine, Bovine, Porcine, Sheep
Conjugation	CF647
Tested Applications	ICC, IF
Concentration	100 µg in 500 µl (0.20 mg/ml)
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Dilution Range	To be determined by the end user
Expiration Date	12 months from date of receipt.