

Product Datasheet

UBXD2 antibody {CF647} orb158679-CF647

Description	Rabbit polyclonal antibody to UBXD2.
Reactivity	Bovine, Equine, Human, Mouse, Porcine, Rat, Sheep
Predicted Reactivity	Canine,Bovine,Porcine
Conjugation	CF647
Tested Applications	ICC, IF
Concentration	100 µg in 500 µl (0.20 mg/ml)
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Dilution Range	To be determined by the end user
Expiration Date	12 months from date of receipt.