

Product Datasheet

ZNF704 antibody {CF647} orb158817-CF647

Description	Rabbit polyclonal antibody to ZNF704.
Reactivity	Human, Mouse, Rat
Predicted Reactivity	Canine, Equine, Porcine
Conjugation	CF647
Tested Applications	ICC, IF
Concentration	100 µg in 500 µl (0.20 mg/ml)
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Dilution Range	To be determined by the end user
Expiration Date	12 months from date of receipt.