

## Product Datasheet

### Rabbit Anti-alpha smooth muscle Actin antibody orb155585-CF647

|                            |                                                                                 |
|----------------------------|---------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to alpha smooth muscle Actin                         |
| <b>Reactivity</b>          | Human, Mouse, Rabbit                                                            |
| <b>Conjugation</b>         | Unconjugated                                                                    |
| <b>Custom conjugations</b> | Unconjugated, Biotin, FITC, HRP, CF647, CF488A, CF594, CF405M                   |
| <b>Tested Applications</b> | ELISA, FC, IF, IHC, WB                                                          |
| <b>Storage</b>             | Shipped at 4?. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| <b>Note</b>                | For research use only.                                                          |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                 |