

Product Datasheet

GPR1 antibody {CF647} orb157289-CF647

Description	Rabbit polyclonal antibody to GPR1.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Predicted Reactivity	Canine
Conjugation	CF647
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from the region aa181-260 of human GPR1
Target	GPR1 {CF647}
Form/Appearance	Liquid: PBS, 2 mM Sodium Azide
Concentration	100 µg (in 411 µl): 0.24 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified by Protein A
MW	41 kDa
Uniprot ID	P46091 (Human)
Entrez	2825
Dilution Range	To be determined by the end user
Expiration Date	12 months from date of receipt.