

Product Datasheet

Cyp11B2 antibody {CF647} orb159192-CF647

Description	Rabbit polyclonal antibody to Cyp11B2 is a member of the cytochrome P450 superfamily of enzymes. It is involved in substrate recognition, conversion and catalyzes the conversion of 11- deoxycorticosterone to aldosterone. This is done by corticosterone and 18- hydroxycorticosterone. Genetic variation in CYP11B2. As a result, they are the cause of corticosterone methyloxidase type 1 deficiency.
Reactivity	Human, Mouse, Rat
Conjugation	CF647
Tested Applications	ICC, IF
Concentration	100 µg in 500 µl (0.20 mg/ml)
Storage	Shipped at 4°. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Dilution Range	To be determined by the end user
Expiration Date	12 months from date of receipt.