

Product Datasheet

Vap1 antibody {CF647} orb159200-CF647

Description	Rabbit polyclonal antibody to Vap1 belongs to the copper/topaquinone oxidase family. It is highly expressed in high endothelial venules of peripheral lymph nodes, appendix, lung and small intestine. It is fundamental in adipogenesis, has adhesive properties, and also has functional monoamine oxidase activity.
Reactivity	Human, Mouse, Rat
Conjugation	CF647
Tested Applications	ICC, IF
Concentration	100 µg in 500 µl (0.20 mg/ml)
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Dilution Range	To be determined by the end user
Expiration Date	12 months from date of receipt.