

## Product Datasheet

### Vnn2 antibody {CF647} orb159208-CF647

<b>Description</b>	Rabbit polyclonal antibody to Vnn2, which is a member of the Vanin family that participate in hematopoietic cell trafficking and in oxidative-stress response. VNN2 functions in transendothelial migration of neutrophils and may control beta-2 integrin-mediated cell adhesion, migration and motility of neutrophil. VNN2 is highly expressed in spleen and blood.
<b>Reactivity</b>	Human, Mouse, Rat
<b>Predicted Reactivity</b>	Canine, Equine, Bovine, Sheep
<b>Conjugation</b>	CF647
<b>Tested Applications</b>	ICC, IF
<b>Concentration</b>	100 µg in 500 µl (0.20 mg/ml)
<b>Storage</b>	Shipped at 4°. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only.
<b>Dilution Range</b>	To be determined by the end user
<b>Expiration Date</b>	12 months from date of receipt.