

Product Datasheet

Wisp3 antibody {CF647} orb159209-CF647

Description	Rabbit polyclonal antibody to Wisp3, which is a member of the CCN family. It is expressed in adult kidney and testis and fetal kidney. It appears to be required for normal postnatal skeletal growth and cartilage homeostasis. Genetic variations in WISP3 are the cause of progressive pseudorheumatoid arthropathy of childhood.
Reactivity	Human, Mouse, Rat
Predicted Reactivity	Canine, Equine, Sheep
Conjugation	CF647
Tested Applications	ICC, IF
Concentration	100 µg in 500 µl (0.20 mg/ml)
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Dilution Range	To be determined by the end user
Expiration Date	12 months from date of receipt.