

Product Datasheet

Ecat1 antibody {CF647} orb159223-CF647

| | |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description | Rabbit polyclonal antibody to Ecat1, which belongs to the KHDC1 family. It appears to be maximal in germinal vesicle oocytes. It may function as a regulator of genomic imprinting in the oocyte. Mutations in ECAT1 are the cause of hydatidiform mole recurrent type 2 |
| Reactivity | Human |
| Conjugation | CF647 |
| Tested Applications | ICC, IF |
| Concentration | 100 µg in 500 µl (0.20 mg/ml) |
| Storage | Shipped at 4°. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only. |
| Dilution Range | To be determined by the end user |
| Expiration Date | 12 months from date of receipt. |