

Product Datasheet

C9Orf61 antibody {CF647} orb159245-CF647

Description	Rabbit polyclonal antibody to C9Orf61, which belongs to the FAM189 family. It prominently expressed in muscle. It is localized to the membrane, Single-pass type I membrane protein.
Reactivity	Human
Conjugation	CF647
Tested Applications	ICC, IF
Concentration	100 µg in 500 µl (0.20 mg/ml)
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Dilution Range	To be determined by the end user
Expiration Date	12 months from date of receipt.