

## Product Datasheet

### HPV31 E7 antibody {CF647} orb159264-CF647

<b>Description</b>	Rabbit polyclonal antibody to HPV31 E7, which can be classified as either high risk or low risk according to their association with cancer. The high risk HPV E7 proteins bind to pRB with a higher affinity than do the low risk HPV proteins.
<b>Conjugation</b>	CF647
<b>Tested Applications</b>	ICC, IF
<b>Concentration</b>	100 µg in 500 µl (0.20 mg/ml)
<b>Storage</b>	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only.
<b>Dilution Range</b>	To be determined by the end user
<b>Expiration Date</b>	12 months from date of receipt.