

Product Datasheet

WBP2 antibody {CF647} orb159268-CF647

Description	Rabbit polyclonal antibody to WBP2, which is an adaptor protein first identified in a screen for proteins that interact with YAP1. WBP2 binds to the WW domain of Yes-associated protein (YAP), WW domain containing E3 ubiquitin protein ligase 1 (AIP5) and WW domain containing E3 ubiquitin protein ligase 2 (AIP2). It interacts with NEDD4.
Reactivity	Human, Mouse, Rat
Predicted Reactivity	Canine, Equine, Bovine, Porcine, Sheep
Conjugation	CF647
Tested Applications	ICC, IF
Concentration	100 µg in 500 µl (0.20 mg/ml)
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Dilution Range	To be determined by the end user
Expiration Date	12 months from date of receipt.