

Product Datasheet

Tetanus Toxin heavy chain antibody {CF647} orb159269-CF647

Description	Rabbit polyclonal antibody to Tetanus Toxin heavy chain, which acts by inhibiting neurotransmitter release. It binds to peripheral neuronal synapses, is internalized and moves by retrograde transport up the axon into the spinal cord where it can move between postsynaptic and presynaptic neurons.
Conjugation	CF647
Tested Applications	ICC, IF
Concentration	100 µg in 500 µl (0.20 mg/ml)
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Dilution Range	To be determined by the end user
Expiration Date	12 months from date of receipt.