

Product Datasheet

ompC antibody {CF647} orb389661-CF647

Description	Rabbit polyclonal antibody to Outer membrane proteins (OMPs) which performs as a permeability channels for nutrients, toxins and antibiotics. Defects in ompC showed a significantly decreased susceptibility to carbapenems and cefepime.
Conjugation	CF647
Tested Applications	ICC, IF
Concentration	100 µg in 500 µl (0.20 mg/ml)
Storage	Shipped at 4°. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Dilution Range	To be determined by the end user
Expiration Date	12 months from date of receipt.