

Product Datasheet

AIF1 (9A3) antibody {CF647} orb704471-CF647

Description	Rabbit monoclonal antibody against AIF1 (9A3) conjugated to CF647
Reactivity	Human, Mouse, Rat
Conjugation	CF647
Tested Applications	ICC, IF
Concentration	100 µg in 500 µl (0.20 mg/ml)
Storage	Shipped at 4°. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Expiration Date	12 months from date of receipt.