

Product Datasheet

Taurine Microplate Assay Kit orb707400

Description	Detection and Quantification of Taurine Content.
Reactivity	All
Range	0.1 mmol/L - 10 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only.
CAS Number	107-35-7
Sample Types	Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
Expiration Date	Please enquire.
Application Notes	Detection Method: Colorimetric, Sample Types: Tissue extracts, Cell lysate, Cell culture media, Other biological fluids, Platform: Microplate reader 400 nm