

Product Datasheet

Coenzyme Q10 Microplate Assay Kit orb707403

Description	Detection and Quantification of Coenzyme Q10 Content.
Reactivity	All
Range	0.1 mmol/L - 10 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only.
CAS Number	303-98-0
Sample Types	Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
Expiration Date	Please enquire.
Application Notes	Detection Method: Colorimetric, Sample Types: Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids, Platform: Microplate reader 620 nm