



Product Datasheet

CD163 antibody orb182468

Description	CD163 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Predicted Reactivity	Canine, Equine, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human CD163 (1001-1121/1156 aa)
Target	CD163
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	130 kDa
Uniprot ID	Q86VB7
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, ICC=1:100-500, IF=1:100-500, Flow-Cyt=1ug/Test, IHC-F=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat colon tissue using CD163 antibody.



Immunohistochemical staining of mouse intestine tissue using CD163 antibody.