

Product Datasheet

SARS-CoV-2 S1 RBD ELISA Kit orb736882

Description	Novel Coronavirus Spike Glycoprotein RBD (SARS-CoV-2 S1 RBD)
Range	1.25 ng/mL-80 ng/mL
Sensitivity	0.31 ng/mL
Target	Novel Coronavirus Spike Glycoprotein RBD (SARS-CoV-2 S1 RBD)
Note	For research use only.
Sample Types	serum, plasma, swabs
Uniprot ID	P0DTC2
Expiration Date	Please enquire.
Application Notes	Detection Wavelength: 450 nm