

Product Datasheet

Human SARS-CoV-2 S1 RBD IgG Antibody ELISA Kit orb736885

Description	Human SARS-CoV-2 S1 RBD IgG Antibody ELISA Kit
Reactivity	Human
Range	Request Information
Sensitivity	Request Information
Target	SARS-CoV-2 S RBD Ab (IgG)
Note	For research use only.
Sample Types	serum, plasma
Uniprot ID	P0DTC2
Expiration Date	Please enquire.
Application Notes	Detection Wavelength: 450 nm