

## Product Datasheet

### Mouse TNFb ELISA Kit orb437531

<b>Description</b>	Mouse Tumor Necrosis Factor Beta (TNFb) Sandwich ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	15.625-1000pg/mL
<b>Sensitivity</b>	5.7pg/mL
<b>Storage</b>	All reagents should be kept according to the instructions on the packaging upon receipt.
<b>Note</b>	For research use only.
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
<b>Uniprot ID</b>	<b>P09225</b>
<b>NCBI</b>	<b>16992</b>
<b>Entrez</b>	<b>16992</b>
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.