

Product Datasheet

Goat Anti-Rabbit IgG(H+L) antibody {CF647} orb688929-CF47

Description	Goat polyclonal antibody against Goat Anti-Rabbit IgG(H+L) conjugated to CF647
Species/Host	Goat
Reactivity	Rabbit
Conjugation	CF647
Tested Applications	ICC, IF
Immunogen	Rabbit IgG(H+L)
Target	Goat Anti-Rabbit IgG(H+L) {CF647}
Form/Appearance	Liquid: PBS and 2mM Sodium Azide
Concentration	100 µg in 500 µl (0.19 mg/ml)
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	The antibody was affinity-purified from goat antiserum by affinity-chromatography using specific immunogen.
Dilution Range	To be determined by end user
Expiration Date	12 months from date of receipt.