

Product Datasheet

PCSK9 Antibody orb745916

Description	PCSK9 Antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB

Immunogen	E.coli-derived mouse PCSK9 recombinant protein (Position: Q182-H560).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	74 kDa
Uniprot ID	Q80W65
Dilution Range	Western blot, 0.25-0.5 µg/ml, Mouse, Rat Direct ELISA, 0.1-0.5 µg/ml, Mouse
Expiration Date	12 months from date of receipt.

Application Notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
--------------------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

