

Product Datasheet

Antibody Purification Kit Protein A orb746992

Description	Proteus Protein A Mini Purification Starter Kit
Tested Applications	Purification
Target	ANTIBODY PURIFICATION KIT PROTEIN A
Preservatives	Buffers contain 0.09% Sodium Azide (NaN ₃) Resin contains 20% Ethanol. MINI - Spin Columns, Buffers, U/F Spinners
Storage	+4°C DO NOT FREEZE
Note	For research use only.
Expiration Date	Please enquire.