

## Product Datasheet

### NF-H antibody orb756864

Description	Gallus polyclonal antibody to NF-H
Species/Host	Gallus
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Custom conjugations	Unconjugated, Biotin, FITC, HRP, CF647, CF488A, CF594, CF405M
Tested Applications	ICC, IHC-P, WB
Immunogen	Native NF-H purified from bovine spinal cord
Target	NF-H
Form/Appearance	Liquid
Concentration	1.0 mg/ml
Storage	Store at 4°C in the dark. Under these conditions, the antibodies should have a shelf life of at least twelve months, provided they remain sterile. Do not freeze these antibodies unless you want to store them for longer periods of time.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	<a href="#">P12036</a>
Dilution Range	WB: 1:2000-1:5000, IHC-P: 1:2000-1:5000
Expiration Date	12 months from date of receipt.