

Product Datasheet

NF-L antibody orb756865

Description	Gallus polyclonal antibody to NF-L
Species/Host	Gallus
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Custom conjugations	Unconjugated, Biotin, FITC, HRP, CF647, CF488A, CF594, CF405M
Tested Applications	ICC, IHC-P, WB
Immunogen	Synthetic peptides
Target	NF-L
Form/Appearance	Liquid
Concentration	300 ug/ml
Storage	Store at 4°C in the dark. Under these conditions, the antibodies should have a shelf life of at least twelve months, provided they remain sterile. Do not freeze these antibodies unless you want to store them for longer periods of time.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P07196
Dilution Range	WB: 1:1000-1:2000, IHC-P: 1:1000-1:2000
Expiration Date	12 months from date of receipt.