

## Product Datasheet

### IBA1 antibody orb756869

|                     |   |
|---------------------|---|
| Description         | Gallus polyclonal antibody to IBA1  |
| Species/Host        | Gallus  |
| Conjugation         | Unconjugated  |
| Custom conjugations | Unconjugated, Biotin, FITC, HRP, CF647, CF488A, CF594, CF405M   |
| Immunogen           | Purified recombinant human full-length IBA1/AIF1 (Uniprot P55008) expressed in E.Coli   |
| Target              | IBA1  |
| Form/Appearance     | Liquid  |
| Concentration       | 0.2 mg/ml   |
| Storage             | Store at 4°C in the dark. Under these conditions, the antibodies should have a shelf life of at least twelve months, provided they remain sterile. Do not freeze these antibodies unless you want to store them for longer periods of time. |
| Note                | For research use only.  |
| Isotype             | IgG   |
| Clonality           | Polyclonal  |
| Uniprot ID          | <b>P55008</b>   |
| Expiration Date     | 12 months from date of receipt.   |