

Product Datasheet

IBA1 antibody orb756870

Description	Gallus polyclonal antibody to IBA1
Species/Host	Gallus
Conjugation	Unconjugated
Custom conjugations	Unconjugated, Biotin, FITC, HRP, CF647, CF488A, CF594, CF405M
Immunogen	Purified recombinant human full-length IBA1/AIF1 (Uniprot P55008) expressed in E.Coli
Target	IBA1
Form/Appearance	Liquid
Concentration	0.2 mg/ml
Storage	Store at 4°C in the dark. Under these conditions, the antibodies should have a shelf life of at least twelve months, provided they remain sterile. Do not freeze these antibodies unless you want to store them for longer periods of time.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P55008
Expiration Date	12 months from date of receipt.