

Product Datasheet

TREM2 antibody orb756872

Description	Gallus polyclonal antibody to TREM2
Species/Host	Gallus
Conjugation	Unconjugated
Custom conjugations	Unconjugated, Biotin, FITC, HRP, CF647, CF488A, CF594, CF405M
Immunogen	Purified recombinant extracellular domain of human TREM2 (Uniprot Q9NZC2) expressed in E.Coli
Target	TREM2
Form/Appearance	Liquid
Concentration	0.2 mg/ml
Storage	Store at 4°C in the dark. Under these conditions, the antibodies should have a shelf life of at least twelve months, provided they remain sterile. Do not freeze these antibodies unless you want to store them for longer periods of time.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NZC2
Expiration Date	12 months from date of receipt.