

## Product Datasheet

### Uricase Microplate Assay Kit orb759230

<b>Description</b>	Detection and Quantification of Uricase Activity.
<b>Reactivity</b>	All
<b>Range</b>	0.03 mmol/L - 3 mmol/L
<b>Tested Applications</b>	FA
<b>Storage</b>	Shipped and store at 4°C for 6 months.
<b>Note</b>	For research use only.
<b>CAS Number</b>	9002-12-4
<b>Sample Types</b>	Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	Detection Method: Colorimetric, Sample Types: Tissue extracts, Cell lysate, Cell culture media, Other biological fluids, Platform: Microplate reader 550 nm