

Product Datasheet

Progesterone Receptor (phospho-Ser190) antibody {CF405M} orb704511-CF405M

Description	Rabbit polyclonal antibody against Progesterone Receptor (phospho-Ser190) conjugated to CF405M
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	CF405M
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from the region around the phosphorylation site of S rat Progesterone Receptor (phospho-Ser190)
Target	Progesterone Receptor (phospho-Ser190)
Form/Appearance	Liquid: PBS and 2mM Sodium Azide
Concentration	100 µg in 400 µl (0.24 mg/ml).
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified by Protein A
MW	99 kDa
Uniprot ID	P06401 (Human); Q00175 (Mouse); Q63449 (Rat)
Entrez	5241
Dilution Range	To be determined by end user
Expiration Date	12 months from date of receipt.