

Product Datasheet

Estrogen Receptor alpha (phospho-Ser106) antibody {CF405M} orb183443-CF405M

Description	Rabbit polyclonal antibody against phospho-Estrogen Receptor alpha (Ser106) conjugated to CF405M
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Predicted Reactivity	Bovine, Equine, Porcine
Conjugation	CF405M
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived between P(p-S)PL amino acids of human Estrogen Receptor alpha (phospho-Ser106)
Target	Estrogen Receptor alpha (phospho-Ser106)
Form/Appearance	Liquid: PBS and 2mM Sodium Azide
Concentration	100 µg in 400 µl (0.24 mg/ml).
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified by Protein A
MW	66 kDa
Uniprot ID	P03372 (Human); P19785 (Mouse)
Entrez	2099
Dilution Range	To be determined by end user
Expiration Date	12 months from date of receipt.