

## Product Datasheet

### Rabbit IgG (Light Chain Specific) antibody orb763499

|                     |                                                                                                         |
|---------------------|---------------------------------------------------------------------------------------------------------|
| Description         | Mouse monoclonal antibody to Rabbit IgG (Light Chain Specific)                                          |
| Species/Host        | Mouse                                                                                                   |
| Reactivity          | Rat                                                                                                     |
| Conjugation         | Unconjugated                                                                                            |
| Tested Applications | IHC, WB                                                                                                 |
| Immunogen           | Recombinant Protein                                                                                     |
| Concentration       | 1mg/ml                                                                                                  |
| Preservatives       | PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles. |
| Storage             | Store at -20°C. Avoid repeated freeze-thaw cycles.                                                      |
| Note                | For research use only.                                                                                  |
| Isotype             | IgG1                                                                                                    |
| Clonality           | Monoclonal                                                                                              |
| MW                  | 25kDa                                                                                                   |
| Dilution Range      | WB 1:5,000-20,000 IHC 1:200-500                                                                         |
| Expiration Date     | 12 months from date of receipt.                                                                         |

