

Product Datasheet

Rat HMGCS ELISA Kit orb782864

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Hydroxymethylglutaryl Coenzyme A Synthase(HMGCS).
Reactivity	Rat
Range	15.63-1000 pg/mL
Sensitivity	5.7 pg/mL
Tested Applications	ELISA
Note	For research use only.
Sample Types	Tissue homogenates and other biological fluids.
Uniprot ID	P17425
Expiration Date	Please enquire.

