

Product Datasheet

Rat LDLR ELISA Kit orb782865

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Low Density Lipoprotein Receptor(LDLR).
Reactivity	Rat
Range	0.32-20 ng/mL
Sensitivity	0.114 ng/mL
Tested Applications	ELISA
Note	For research use only.
Sample Types	tissue homogenates, cell lysates and other biological fluids
Uniprot ID	P35952
Expiration Date	Please enquire.

