

## Product Datasheet

### Fish MT1 ELISA Kit orb782867

<b>Description</b>	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Metallothionein 1(MT1).
<b>Reactivity</b>	Zebrafish
<b>Range</b>	62.5-4000 pg/mL
<b>Sensitivity</b>	26.3 pg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only.
<b>Sample Types</b>	serum, plasma, tissue homogenates and other biological fluids
<b>Uniprot ID</b>	<b>P52722</b>
<b>Expiration Date</b>	Please enquire.

