

Product Datasheet

Rat Cabin1 ELISA Kit orb782868

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Cabin1(Calcineurin-binding protein cabin-1).
Reactivity	Rat
Range	62.5-4000 pg/mL
Sensitivity	26.3 pg/mL
Tested Applications	ELISA
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates and other biological fluids
Uniprot ID	O88480
Expiration Date	Please enquire.

