

## Product Datasheet

### CA125 antibody {HRP} orb10214-HRP

Description	Rabbit polyclonal antibody against CA125 conjugated to HRP
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	HRP
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human CA125. Please <a href="#">contact us</a> for the exact immunogen sequence. The peptide is available as <a href="#">orb374704</a> .
Target	CA125
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	100 µg (in 522 µl): 0.192 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	2353 kDa
Uniprot ID	<a href="#">Q8WXI7</a>
Dilution Range	To be determined by the end user
Expiration Date	12 months from date of receipt.