

Product Datasheet

IL6 antibody {CF647} orb665927-CF647

Description	Rabbit polyclonal antibody against IL6 conjugated to CF647
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	CF647
Tested Applications	ICC, IF
Immunogen	Recombinant protein of human IL-6. The exact sequence is proprietary.
Target	Recognizes endogenous levels of IL-6 protein.
Form/Appearance	Liquid: PBS and 2mM Sodium Azide
Concentration	100 µg in 400 µl (0.24 mg/ml)
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	The antibody was purified by immunogen affinity chromatography.
Uniprot ID	P05231 (Human); P08505 (Mouse); P20607 (Rat)
Entrez	3569
Dilution Range	To be determined by end user
Expiration Date	12 months from date of receipt.