

Product Datasheet

Cadherin-17 antibody {HRP} orb11244-HRP

Description	Rabbit polyclonal antibody against Cadherin-17 conjugated to HRP
Reactivity	Human, Mouse, Rat
Conjugation	HRP
Tested Applications	ICC, IF
Concentration	100 µg in 500 µl (0.20 mg/ml)
Storage	Shipped at 4°. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Dilution Range	To be determined by the end user
Expiration Date	12 months from date of receipt.