

Product Datasheet

RAB5B antibody {Biotin} orb652383-Biotin

Description	Mouse monoclonal antibody against RAB5B conjugated to Biotin
Species/Host	Mouse
Reactivity	Human, Monkey
Conjugation	Biotin
Tested Applications	ICC, IF
Immunogen	This RAB5B antibody is generated from a mouse immunized with a recombinant protein of human RAB5B.
Target	RAB5B
Form/Appearance	Liquid: PBS and 2mM Sodium Azide
Concentration	100 µg in 400 µl (0.24 mg/ml)
Storage	Aliquot and store at -20°C or -80°C. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG1,k
Clonality	Monoclonal
Purity	Affinity purified
MW	23707 Da
Clone ID	1615CT175.80.21
Uniprot ID	P61020
Entrez	5869
Dilution Range	To be determined by end user
Expiration Date	12 months from date of receipt.