



Product Datasheet

TIMP-1 antibody orb313247

Description	TIMP-1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rabbit
Predicted Reactivity	Bovine, Canine, Porcine, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human TIMP-1 (103-207/207 aa)
Target	TIMP-1
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	21 kDa
Uniprot ID	P01033
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human colon carcinoma tissue using TIMP-1 antibody.



Immunohistochemical staining of human gastric carcinoma tissue using TIMP-1 antibody.